

We advise you to store this Bulletin with the applicable User Manual.

Smart View or SmartClinic may not Trigger Calculations for BSA, CrCl or GFR

Product: Smart View or SmartClinic and MOSAIQ®

Date: October 2024

Reference number (Field Change Order, FCO): 371-05-MSQ-081

Field Corrective Action number (FCA), if applicable: NA

Relates To: SVJ-11230, SVJ-11470, SVJ-11471

This document contains confidential and proprietary information of the Elekta Group and is for recipient use only. Subject to Copyright protection, any dissemination, distribution or copying of this document is strictly prohibited without written permission of Elekta.

SVJ-11230:

Scope:	MOSAIQ 3.0 and above using Smart View or SmartClinic.
Description:	When an assessment includes GFR values and the patient's height, weight, creatinine, BUN, and Cystatin C are entered in Smart View or SmartClinic, clicking Save may trigger an error, preventing the calculation of GFR(Schwartz), GFR(CKD-EPI), and GFR(Cystatin). In this case, MOSAIQ will revert to the last calculated GFR values.
Solution:	After acknowledging the error message and confirming the GFR is not calculated in Smart View or SmartClinic, enter the input values for GFR calculations directly into MOSAIQ. This defect will be resolved in a future Smart View release.
Technical reference:	NA
Contact:	If you have any queries about this Notice, please contact your local Elekta office.

SVJ - 11470

Scope:	MOSAIQ using Smart View or SmartClinic.
Description:	If an assessment view is set up to include height and weight but omits BSA, entering height and weight in Smart View or SmartClinic will not trigger a BSA calculation. MOSAIQ will use the last known calculated BSA value.
Solution:	If entering height and weight in Smart View or SmartClinic, ensure the assessment view contains BSA. This defect will be resolved in a future Smart View release.
Technical reference:	NA
Contact:	If you have any queries about this Notice, please contact your local Elekta office.

SVJ - 11471

Scope:	MOSAIQ using Smart View or SmartClinic.
Description:	When using Smart View or SmartClinic in a translated environment, users may encounter a “missing gender” error for CrCl and GFR values that require gender for their calculations, despite the patient having a gender recorded. In such cases, MOSAIQ will default to previously calculated CrCl and GFR values.
Solution:	After acknowledging the error message and confirming the CrCl or GFR are not calculated in Smart View or SmartClinic, enter the input values for CrCl or GFR calculations directly into MOSAIQ. This defect will be resolved in a future Smart View release.
Technical reference:	NA
Contact:	If you have any queries about this Notice, please contact your local Elekta office.